## Finding your serial number on a LOGIQ System

On touch panel select Utility> Admin

System Admin Terrs Logon Groups System Password	Disk Encryption	Audit
roduct	Option	Sti
Product Radiology Eagle	Co	odedContrast Perma
HW Number 503061US7		
SW Option Key	Shear Wave	Elastography Permanent
Enter New Option Key Add		CW Permanent
	Automated Funct	ional Imaging Permanent
Installed Option Keys 6YT87-APMZ7-NPM34-KV2TM-JDGKA		
UCRQE-WJZHL-%NENS-H8ANH-26ARL		BSteer+ Permanent
Remove		AutoIMT Permanent
Import		
Service		
Enable Automatic Request for Service 🥃		
Transfer Machine Logs to GE 🛛 😨 Transfer Operator Information to GE		
Enable Automatic Probe Diagnostics		
Transfer Images (captured without patient information) to GE		
Include Annotations with Image Transfer to GE		
Protected Health Information(PHI)		
Prevent Writing to Removable Media-CD/DVD/USB (requires reboot)		
USB External Media		
USB External Media disabled (requires reboot)		
By Checking this box you will disconnect all external USB Mass Storage devices.		
Rights		
Require Admin Operator Rights to Save Imaging Settings 🔲		
Save Exit Search Cancel		